



CARE HEALTH CARE LIFE
Dr.@home Series

LIGHT THE RHYTHM OF LIFE
WITH TECHNOLOGY



Portable Blood Pressure Monitor DX10

V6CO. Limited Liability Co.

1185 Avenue of the Americas 3rd Floor
New York NY 10036 United States

1-833-838-V6CO (8626)

support@v6co.com

www.v6co.com

V6CO. reserves the right of product appearance, specifications and function modification and updating, without prior notification.

V6CO. Copyright 2022 V6CO. Limited Liability Company. All rights reserved.

Portable Blood Pressure Monitor

DX10



DX10 is an LED screen Electronic Sphygmomanometer, simple interface and one-key measurement, it can measure NIBP accurately, NIBP records can be reviewed by the memory button and support voice broadcasting(optional). The device is widely used in family, clinic and physical-examination center for routine test.

Performance

- ◆ Measurement method: Oscillometry
- ◆ Measurement mode: upper-arm type
- ◆ Measurement range: 0 kPa(0 mmHg) ~ 39.6kPa(297 mmHg)
- ◆ Resolution: 0.133kPa(1 mmHg)
- ◆ Accuracy: +0.4kPa(+3 mmhg)
- ◆ PR measurement range: 40 bpm ~ 240 bpm
- ◆ Inflation: automatic Inflation by force pump
- ◆ Deflation: automatic multistep deflation
- ◆ Display: LED large colour screen
- ◆ Power: Four "AA" batteries/5V power adapter
- ◆ Safety character: Class II Type BF

Physical Characteristic

Dimension:130mm(L)x110mm(W)x80mm(H)

Operation environment:

Temperature:+5c~+40c

Relative humidity: 15%~85%

Atmospheric pressure: 700 hPa~1060 hPa

Storage environment:

Temperature:20°C ~ 55°C

Relative humidity: < 95%

Atmospheric pressure: 700 hPa ~1060hPa

Features

- ◆ Small volume, one- key measurement, easy to operate, big font in displaying, simple and easy to understand the displayed contents.
- ◆ Support voice broadcasting (optional).
- ◆ Start measuring manually, record each measurement data and store up to 199 groups of data.
- ◆ LED colour screen. Automatic power off when there is no operation for a long time, achieves saving power.
- ◆ High pressure prompt.
- ◆ Low battery and error information indication.
- ◆ Optional units: mmHg and kPa.

Accessories

- ◆ Cuff for adult
- ◆ User manual

