




CARE HEALTH CARE LIFE
Dr.@home Series



DESIGNED FOR CHILDREN



EDUTAINMENT



DO IT YOURSELF



SAFE



KT20

Kid Pulse Oximeter

V6CO. Limited Liability Co.

1185 Avenue of the Americas 3rd Floor
New York NY 10036 United States

1-833-838-V6CO (8626)

support@v6co.com

www.v6co.com

V6CO. reserves the right of product appearance, specifications and function modification and updating, without prior notification.
V6CO. Copyright 2022 V6CO. Limited Liability Company. All rights reserved.




KT20

Kid Pulse Oximeter




Technical Parameters

- ◆ Integrated with SpO₂ probe and processing display module
- ◆ Application: hospital, home, community medical treatment, sports healthcare, etc
- ◆ Can measure SpO₂ and Pulse Rate accurately
- ◆ SpO₂ and Pulse Rate display, Pulse Rate waveform and bar graph display
- ◆ Battery voltage low indication
- ◆ The display mode can be changed
- ◆ Screen brightness adjustable
- ◆ Low power consumption
- ◆ Small volume(59(L) × 36(W) × 35(H) mm)、light, easy to carry
- ◆ Display Mode:double Color display
- ◆ SpO₂ Measuring Range:0 ~ 100%
- ◆ PR Measuring Range:30 ~ 250bpm
- ◆ Resolution: 1% for SpO₂, 1bpm for pulse rate
- ◆ Accuracy: ±2% (70%~ 100%), unspecified (<70%)
- ◆ for SpO₂, ±2bpm or ± 2% (select larger) for pulse rate
- ◆ Power Supply: 1.5V (AAA size) alkaline batteries×2
- ◆ Can work 36 hours continuously
- ◆ Interference resistance capacity against ambient
- ◆ Measurement performance at low perfusion



With DIY stickers, children can collage as they like, making traditional medical products more interesting.



CMS50Q1

Enjoyment

Innovative children's medical testing product experience

