



**INNOVATIVE,
TECHNOLOGY
& WELLNESS**



V6CO was founded in 2019 and began to provide wellness devices to people worldwide. The company consists of a small and innovative team of creators, engineers, designers, and dreamers passionate about improving the quality of life for people with medical conditions that require long-term care at an affordable price. V6CO is based in the USA and has a growing number of customers worldwide.

The COVID-19 pandemic changed how people seek medical care, with more preferring in-home services and V6CO is devoted to meeting these health care demands. Its wellness devices are developed and produced with the latest technologies. The medical-grade components work together seamlessly, meaning less time at the doctor's or in long-term nursing or hospital care.

The home medical device startup's mission is to provide quality long-term care products that improve the quality of life for those who need it most. More people are living longer lives through the advancement of medical technology. V6CO's goal is to promote wellness for these people who require long-term medical care at home.

From hearing devices to mobility aids, their products are designed to help people with chronic conditions live fuller, healthier lives in the comfort of their own homes. Customers can purchase the medical-grade, life-changing technology at their local CVS or Costco.

HI! BABY SOUNDER FETAL DOPPLER



PORTABLE FETAL DOPPLER HI! BABY SOUNDER.

Proudly cleared by FDA 510(k), with
precise accuracy testing result.



PERFORMANCE

PORTABLE FETAL DOPPLER

Safety classification:

internally powered equipment, type
CF applied part

Display:

1.77 " TFT screen

FHR measurement and display range:

50 BPM ~ 240 BPM (BPM: beats per minute)

Resolution: 1 BPM**Accuracy:** ± 2 BPM**Rated working voltage:** DC3 V**Effective Radiating Area of Transducer:**

$\leq 208 \text{ mm}^2$

Battery type: two 1.5V batteries (AA LR6)**Nominal Frequency:** 3.0 MHz**Working Frequency:** 3.0 MHz $\pm 10\%$ **Negative Peak Sound Pressure:** $p_- < 1 \text{ MPa}$ **Output Beam Intensity:** $I_{ob} < 20 \text{ mW/cm}^2$ **Spatial-peak Temporal-average Derived Intensity:**

$I_{spta} < 100 \text{ mW/cm}^2$

Ultrasonic Output Power: $P < 20 \text{ mW}$ **Working Mode:**

Continuous wave Doppler





FEATURE & BENEFITS

◆ LIGHT WEIGHT & EASY TO HOLD

◆ CURVED WAND DESIGN

◆ PRECISION ENGINEERING

◆ HEAR YOUR BABY

◆ WATERPROOF

◆ CRYSTAL CLEAR SCREEN

DR.@HOME DOPPLE GELIXIR



| GELIXIR

Mega Hydrating

Hyaluronic Acid & Vitamin E are powerful natural humectants that provide long-lasting hydration to improve the look of skin texture.

Calming & Repairing

Arnica Montana Flower contains natural lactone compounds & helenalin. Anti-inflammatory, Antimicrobial, Astringent & Antiseptic. Help to reduce redness & itchiness.

Anti-aging

Powerful α -Tocopherol Acetate helps skin to self-repair damage, and improve the appearance of spots, discoloration and signs of aging.

Coupling Agent

FDA Approved Medical Grade Ultrasound Gel with Low Allergic Rate (Different Viscosity)





What are the ingredients

natural ingredient

in your Dopple Gelixir?



α -Tocopherol
Acetate



Betula Alba (Birch)
Bark Extract



Carbus Arnica
Extract



Centella Asiatica
Extract



Calendula
Officinalis Extract



Concentrated
Hyaluronic Acid

5-IN-1 MULTIFUNCTION SEATED CANE



5-IN-1 WALKING CANE WITH FOLDING SEAT



Lightening & USB

LED light
SOS light
USB charging station

FM Radio

Easy access control panel

SOS Alarm System

When accidentally fall,
alarm system will detect and
attract people around to help

Joint Support Seat

Rivets and screws fixed, seat plate and
aluminum rod connection more firmly

Height Adjustable

Adjustable Height, adapt to people of
different heights





5-IN-1 WALKING CANE WITH FOLDING SEAT

SGS test result - seat loading capacity ~ 440lbs,
handle capacity ~ 330lbs

Hard to find similar level products on the market that can take more than 200lbs with stable seated feature. Plus, extra 5 addition functions bring lots of competition power to the product, yet price lower than majority of the lower quality product on the market now.



GUARDIAN AIR VEST



GUARDIAN AIR VEST


The airbags deploy to protect the body with algorithms recognize changes in user's states within 0.08s.

Using inserted gas cylinder, compared with carbon dioxide gas that commonly used in the market, Helium Gas is much more stable and high temperature resistance. Our high quality Helium Gas Cylinder can tolerate temperature range from lowest -60°C to highest 85°C .



GUARDIAN AIR VEST FEATURES

TECHNICAL INFORMATION



| FEATURES

Inflator; High Precision Microchip Sensor

Charging Current: 2A

Charging Voltage: 110-230V

Large Battery Life Capacity: 5000 mAh

(Fully charged can last 8 days)

INCLUDES

Guardian Air Vest - X1

USB Power Adapter

USB Magnetic Charging Cable

User Manual



KID PULSE OXIMETER





| KIDS PLUSE OXIMETER

Proudly cleared by FDA 510(k), a medical grade accuracy device

For kids 3-12 years old. Aim at children's size, requirements and preference. Easily wear on one finger. Convenience for daily use. Record oxygen levels & heart rate every second continuously



**ONE BOX OF HEALTH,
TWO TYPES OF FUN**

**KIDS PULSE
OXIMETER**



KID PLUSE OXIMETER



DESIGNED
FOR CHILDREN



EDUTAINMENT



DO IT YOURSELF



SAFE





PORTABLE BLOOD PRESSURE MONITOR

PORTABLE BLOOD PRESSURE MONITOR

FEATURES

Small volume, one-key measurement, easy to operate, big font in displaying, simple and easy to understand the displayed contents.

Support voice broadcasting (optional).

High pressure prompt.



Start measuring manually, record each measurement data and store up to 199 groups of data.

LED colour screen. Automatic power off when there is no operation for a long time, achieves saving power.

Low battery and error information indication.



PORTABLE BLOOD PRESSURE MONITOR

DX10 is an LED screen Electronic Sphygmomanometer, simple interface and one-key measurement, it can measure NIBP accurately, NIBP records can be reviewed by the memory button and support voice broadcasting(optional). The device is widely used in family, clinic and physical-examination center for routine test.





SMART ECG MONITOR


SMART ECG MONITOR

Portable ECG Monitor is a device for checking ECG, which is applicable for family and individual user, it is a good helper to early prevent from cardiovascular diseases and reduce risks.

Intelligent design, achieves remote health management by using with mobile application, it can automatically start measurement, store ECG, upload data, download health conclusion and obtain doctor advice at the first moment.




PORTABLE ECG MONITOR



PERFORMANCE

- ◆ Calibration voltage : $1\text{mV} \pm 5\%$
- ◆ Standard sensitivity : $10\text{mm/mV} \pm 5\%$
- ◆ Amplitude frequency characteristic : standard: $10\text{Hz}, 1\text{Hz} \sim 20\text{Hz}$; $(+0.4\text{dB}, -3\text{dB})$
- ◆ Noise level : $\leq 30\mu\text{V}$
- ◆ CMRR : $\geq 60\text{dB}$
- ◆ Scanning speed : $25\text{mm/s} \pm 5\%$
- ◆ Sampling rate : 250 dots/s
- ◆ HR measurement range : $30\text{bpm} \sim 300\text{bpm}$, error : $\pm 1\text{bpm}$ or 1%
- ◆ Type of protection against electric shock: Internal power device
- ◆ Degree of protection against electric shock: Type BF applied part
- ◆ Degree of waterproof: IP22





AIR SOUNDER HEARING AIDS

AIR SOUNDER HEARING AIDS

Hear the conversations you want without the hassle of going through a doctor or audiologist for a fitting.

Air Sounder Hearing Aids provides a full-digital smart hearing solution without paying for a doctor's exam or fitting costs. The user can carry out a listening test independently via APP and then generate the optimal hearing compensation scheme according to their own needs and preferences.



COMPARE V6CO AIR SOUNDER TO OTHER

	V6CO. AirSounder	BOSE SoundControl
Peak Output OSPL90	114 dB	113 dB
HF Average OSPL90	106 dB	106 dB
Peak Full-on Gain	36 dB	N/A
HF Average Full-on Gain	30 dB	30
WDRC Frequency Range	200 ~ 6000Hz	200 ~ 8000Hz
Total Harmonic Distortion 500Hz	1.10%	1%
Total Harmonic Distortion 800Hz	0.40%	0.50%
Total Harmonic Distortion 1600Hz	0.40%	0.50%
Equivalent Input Noise Level	20 dB	27 dB
Battery Current	5mA	2.8mA
Price Range	\$\$	\$\$\$
	\$668USD	\$849USD

AIR SOUNDER HEARING AIDS

Set Using APP Listening Test

Skip the doctor's exam and professional fitting. Air Sounder hearing aids enable you to manage your own personal sound profile.

A Multifunction Design

Take all your calls with the in-ear microphone and listen to music without switching devices.




Affordable

The Air Sounder hearing aids provide the best quality at a reasonable price.

Simple and Smart Design

Air Sounder uses smart algorithms to monitor your listening environments and adjust the sound perfectly for you, so you can hear better wherever you go.

TINNITUS CORRECTION AIDS



TINNITUS CORRECTION AIDS

Regulated by FDA 510(k) Hospital Grade

Dr.@home series tinnitus hearing aids, skip a medical test, directly to wear" smart digital hearing aids, utilizing the latest US nano chip digital tech, are priced competitively, top notch performance



HEALTH MOBIOUS MONITORING RING



HEALTH MONITORING RING

Epoxy Crystal Ceramic, Constantly Monitoring
Body Temperature, BPM
Constantly Monitoring, Fast Access, etc.



HEALTH MONITORING RING



DIABETES SOCKS



DIABETES GRIPPED SOCKS

Keep Your Feet Safe and Comfortable

Diabetes Grippped Socks are designed to protect your feet and keep you safe while living with diabetes. The socks have a fully padded under-sole for active or sports styles, as well as treads to prevent falls



DIABETES GRIPPED SOCKS

Breathable Material

The tight grip crew sock is made of Circufiber, which is breathable and great for sensitive skin.

Hand-Sewn Heel Seam

The precise hand-sewn heel seam is free from irritation and there is no compression applied.

Slip-Free Grip

The unique grips on the bottom of the sock keep you from slipping around and firmly on your feet.

Odor-Free Material

The sock features Dri-Release Geo, which is moisture-wicking and anti-microbial for odor control.

Non-Binding Top Band

The top band on the sock is non-binding and will not create any irritation or compression.

Fully Padded Under Sole

The sock has a fully padded under-sole for added protection and support as well as adding warmth.





PORTABLE AIR PURIFIER NECKLACE

PORTABLE AIR PURIFIER NECKLACE

Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua. Ut enim ad minim veniam.



PURIFICATION EXPERIMENT EFFECT



01 Prepare glass, cigarette, lighter, air purifier and other materials



02 Light the cigarette to fill the glass with smoke, and put the starting air purifier in the left cup



03 Observe the smoke changes of the tow cups. After a few minutes, the smoke of the left cup is completely decomposed.



**MEANINGFUL INNOVATIONS
FOR EVERYDAY LIFE,
AVAILABLE WHEN & WHERE YOU NEED IT.**

